

Step 1:



Draw 10mL of Sodium Citrate Anticoagulant into 60mL Syringe

Step 2:



Draw 50mL whole blood from the patient, filling the syringe to 60mL

Step 3:

Please remove and **DISCARD RED CAP** from concentrating device before use



Load anticoagulated whole blood into the **Concentrating Device**


Step 4:



Counterbalance and process the **Concentrating Device** at

5 minutes
4400 RPM


Step 5:



Piston touches the RBC interface
(trapping plasma & platelet buffycoat inside)

Attach the 60mL syringe and aspirate the plasma until the piston touches the RBC interface (trapping the plasma & platelet buffycoat inside), then stop aspirating

Step 6:



7mL PRP
(Plasma & Platelet buffycoat)

Attach the 12mL syringe and aspirate 7mL PRP (Rotate the syringe to re-suspend the platelet buffycoat into the plasma)